

This tool is intended to be used by staff/consultants installing GxAlert remote monitoring system in a GeneXpert equipment.

The Installation Checklist is used to ensure correct installation of the GxAlert system. All sections must be filled and then scanned and send to GeneXpert focal person. Make sure not tests are running on GeneXpert. If they are, wait until they complete before starting install on computer (contact lab in advance).

1. Information						
Date of installation				Name of installer:		
Facility name/ Laboratory name and GX serial number						
2. Install NTLP Secure Router in secure location						
Router serial number				Where was it installed?		
3. Power up secure router and connect to wireless network. Check active SIM and tick						
Active SIM card	SIM 1 <input type="checkbox"/>	2G <input type="checkbox"/>	Signal strength ____%	Verify if it says VPN CONNECTED		YES <input type="checkbox"/>
	SIM 2 <input type="checkbox"/>	3G <input type="checkbox"/>				NO <input type="checkbox"/>
Install wireless adapter		YES <input type="checkbox"/>	Access https://admin.gxalert.com	YES <input type="checkbox"/>	If no, check and repeat all	
		NO <input type="checkbox"/>		NO <input type="checkbox"/>		
4. Start GeneXpert Software						
Configure Host Communication		YES <input type="checkbox"/>	Review completed and Upload expired tests for the current month:		YES <input type="checkbox"/>	
		NO <input type="checkbox"/>			NO <input type="checkbox"/>	
Set system to Automatic Result Upload:		YES <input type="checkbox"/>	Install GXConnect:		Verify SMS or Email that site is connected	
		NO <input type="checkbox"/>				
					NO <input type="checkbox"/>	
5. Record Cartridge supply at facility by expiration date						
Number of cartridges:				Exp date		
Number of cartridges:				Exp date		
Number of cartridges:				Exp date		